

Nu Seal

FLUORIDE RELEASING LIGHT-CURED CAVITY LINER



Nu Seal has a flowable consistency that facilitates reinforcing cavity floors and lining cavity walls prior to application of filling materials. This material consists of a biocompatible light-cured resin, filled with hydroxyapatite, barium sulfate and fluoride salts.

Features:

- Fast and easy to apply, convenient direct-to-tooth application.
- Mechanically strong (over 15,000 psi, 104 MPa), virtually insoluble in oral fluids.
- Fluoride releasing to strengthen the dentin and prevent decay.
- Compatible with all commercially available restoratives and cements.
- X-ray opaque.

Four Syringe Kit: Item # 914
4 - 1.5 gm syringes, 20 prebent applicator tips and instructions.

One Syringe Kit: Item # 907
1 - 1.5 gm syringe, 5 prebent applicator tips and instructions.

Additional Tips: Item # 012
100 - 20 gauge prebent applicator tips.

*Because patients will continue
to whiten for a lifetime!*

Nu Radiance

Forté[®]
with
CaCP

**Advanced Teeth Whitening
that restores enamel health and
improves luster!**



Nu Radiance Restoratives are available through most dental distributors.

For more information, please visit our website at:
www.nuradiance.com

Nu Radiance, Inc.
Green Bay, WI 54304
Toll Free 866-899-3207
Tel. 920-965-0953
Fax. 920-965-0954
Email: info@nuradiance.com

Nu Radiance
RESTORATIVES



Nu Form™

Temporary Crown & Bridge Material 1:1
With Excellent Flexural Characteristics



Nu Form is a composite based, self-curing resin material designed for easy fabrication of temporary crowns and bridges. It sets in two stages, providing for easy working and exceptional durability. For making repairs, Nu Form is compatible with self and light-cure resin based material.

Features:

- Initial set 70 seconds!
- Final set 4.0 minutes!
- Phenomenally tough.
- High resistances to impact.
- Fluorescence assures good appearance in all lighting conditions.
- Minimal oxygen inhibited layer.
- Easy to polish.
- Transient flexibility, provides ease of handling.

Kit Contains: 50 ml of material in an auto-mix cartridge 1:1, 15 mixing tips and instructions.

Shades Available:

Bleach	Item # 808	A1/B1	Item # 709
A1	Item # 761	A2	Item # 716
A3	Item # 723	C2	Item # 730

Additional Tips: Item # 747
50 blue mixing tips per bag.

Nu Tempcement NE

Non-Eugenol Polymeric Temporary Cement



Nu Tempcement NE is a two phase resin based temporary dental cement and luting agent. Nu Tempcement cements temporary restoratives (temporary crowns and bridges) while the permanent prosthesis is being made and to assess if any trauma resulting from abutment's preparation has subsided.

Features:

- Excellent handling characteristics.
- X-ray opaque.
- Negligible solubility.
- Low film thickness.
- Compatible with all commercially available permanent cements.
- Well-balanced adhesive properties provide good prosthesis retention.

Two Syringe Kit: Item # 778

2 - 5.5 gm hand held auto-mix syringes of Nu Tempcement, 15 auto-mix tips, 1 - 3 gm syringe of Release Agent and instructions.

One Syringe Kit: Item # 785

1 - 5.5 gm hand held auto-mix syringe of Nu Tempcement, 8 auto mix tips, 1 - 3 gm syringe of Release Agent and instructions.

Additional Tips: Item # 792
50 brown mixing tips per bag.



1st Core is a self-curing, auto-mix core build-up material. This smooth flowing, non-slumping formula offers great handling characteristics; consistent working time and fast, convenient, full-depth curing.

Features:

- Great handling characteristics and consistent working time (minimum 60 sec.)
- Fast (120 seconds), full-depth curing - Reduces chair time!
- Great mechanical strength (ca 36,000 psi.)
- Excellent x-ray opacity.
- Available in three desirable shades.

Kit Contains: 68 gm (~48 ml) of material in an auto-mix cartridge 1:1 ratio, 15 mixing tips, 15 applicator tips and instructions.

Shades Available:

White	Item # 419
Natural	Item # 402
Contrast	Item # 426

Additional Tips: Item # 400
50 mixing tips and 50 applicator tips per bag.